

## SASHIMI hotkeys and mouse actions

Keyboard	Action
Arrow left/right	Scroll down/up slices.
d	Delete a point – click on a point to delete it.
a	Add object – click on the image to create a polygon object.
r	Remove object – click inside object to remove it.
t	Track objects.
u	Update last tracking task with updated settings.
m	Toggle mask visibility.
s	Toggle object (polygon+spline) visibility.
. (dot)	Increase mask opacity.
, (comma)	Decrease mask opacity.
l	Toggle label name visibility.
1-9	Select label 1-9.
Ctrl + o	Open image.
Ctrl + l	Open object file.
Ctrl + s	Save object file.
Ctrl + e	Export segmentation mask as NIfTI file.
Ctrl + z	Undo last action.
Ctrl + y	Redo last action.

Mouse	Action
Scroll wheel	Zoom in/out.
Ctrl + scroll wheel	Scroll down/up slices.
Hold scroll wheel + move	Pan axes.